



## **RCA120, AGAROSE BOUND**

**SKU:** AL-1083-5



---

### **DESCRIPTION**

Agarose bound *Ricinus communis* agglutinin I is prepared using our affinity-purified lectins. This lectin consists of two subunits of 60 kDa which can be dissociated by reducing agents into closely related chains between 27 kDa and 33 kDa. One of the chains appears to be common to the B chain of another castor bean lectin, ricin, while the other chain is unique to RCA I.

#### **Features:**

- Bead diameter ranges in size from 45-165 microns
- Matrix is stable in solutions at pH 3-11 as well as many organic solvents
- Immobilized lectins are prepared using affinity purified lectins
- Covalent attachment preserves lectin activity and minimizes conformational changes that might result in nonspecific or hydrophobic interactions
- Conjugated proteins are not leached off the beads by Tris or other routinely used buffers
- No residual charges present after conjugation. This minimizes non-specific binding to the matrix

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



- Product supplied as a 1:1 suspension in buffer
- Inhibiting/Eluting Sugar: 200 mM galactose or lactose or Glycoprotein Eluting Solution (ES-2100)

## SPECIFICATIONS

<b>Molecular Weight</b>	120
<b>Extinction Coefficient</b>	1.17
<b>Formulation</b>	10 mM HEPES, pH 7.5, 0.15 M NaCl, 20 mM lactose, 0.08% sodium azide
<b>Inhibiting or Eluting Sugar</b>	Lactose
<b>Label Modifier Type</b>	Lectins
<b>Unit Size</b>	5 ml
<b>Storage Instructions</b>	2-8 °C DO NOT FREEZE
<b>Sugar Specificity</b>	Terminal type 2 LacNAc
<b>Usage Summary</b>	Wash gel thoroughly with buffer before use to remove sugar added to stabilize the lectin. Recommended product for eluting glycoconjugates bound to this agarose-lectin: Glycoprotein Eluting Solution, Cat. No. ES-2100. Alternatively, 200 mM galactose or lactose can be used. After use, wash the gel with several column volumes of buffered saline, then resuspend gel in buffered saline containing 0.08% sodium azide for storage.
<b>Applications</b>	Glycobiology, Affinity Chromatography
<b>Conjugate</b>	Agarose

## TECHNICAL INFORMATION

Agarose bound\* *Ricinus communis* agglutinin I is prepared using our affinity-purified lectins. Heat stable, cross-linked 4% agarose beads with a molecular weight exclusion limit of about  $2 \times 10^7$  daltons are used as the solid-phase matrix to which the lectins are covalently coupled. The attachment of the lectins to the beads is carefully controlled to preserve lectin activity and minimize conformational changes of the bound lectins that might result in nonspecific ionic or hydrophobic interactions. The technique we have developed to couple lectins to agarose beads

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



inserts a hydrophilic spacer arm between the lectin and the matrix.

This coupling method provides several advantages over the traditional cyanogen bromide procedure:

- Maximum carbohydrate binding activity of the coupled lectins is retained
- Linkage is stable over a range of pH values

Our agarose bound lectins are supplied at a constant concentration of lectin per ml of settled beads. The concentration for each lectin is selected to achieve the highest glycoconjugate binding capacity per mg of lectin present in the beads. Each lot is tested for its binding capacity using glycoproteins known to bind the lectin. This provides a guideline for the user and assures the quality of our agarose bound lectins.

\*4 mg lectin/ml gel

## DOCUMENTS

- [Safety Data Sheet](#)
- [Lectins in Histochemistry, ELISA, and Western Blot Applications](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**