


AZIDO-DPEG®₈-NHS ESTER

SKU: QBD-10503



Azido-dPEG®₈-NHS ester, product number QBD-10503, is a water-soluble click chemistry crosslinker containing an azide group linked to an N-hydroxysuccinimidyl (NHS) ester through a single molecular weight, discrete polyethylene glycol (dPEG®) spacer. Azido-dPEG®₈-NHS ester works with copper(I)-catalyzed, ruthenium-catalyzed and with copper-free (e.g., strain-promoted) click chemistry. The dPEG® spacer imparts water solubility and adds hydrodynamic volume to the conjugated product. The single molecular weight product design, with its discrete chain length, simplifies the analysis of this product.

Our single molecular weight, discrete chain length PEG (dPEG®) products are a superior alternative to traditional polyethylene glycol products. Conventional PEGylation reagents consist of a complex mixture of different chain lengths due to the polymerization processes that form them. In contrast, each dPEG® product contains a discrete chain length of PEG with only one molecular weight. Conjugates made with dPEG® products have less complicated analyses compared to conjugates prepared with traditional PEGylation reagents.

NHS esters are the most popular, most widely used way to conjugate carboxylic acids to primary or secondary amines resulting in stable amide bonds. NHS esters react quickly and efficiently in aqueous media at physiological pH values (7.0 – 7.5). However, they are prone to hydrolysis at a rate that is pH-dependent. Published research shows that 2,3,5,6-tetrafluorophenyl (TFP) esters are more hydrolytically stable and have better reactivity with amines than NHS esters.

Specifications

Unit Size	100 mg, 1000 mg
Molecular Weight	564.58; single compound
Chemical formula	C ₂₃ H ₄₀ N ₄ O ₁₂

For research use only. Not intended for animal or human therapeutic or diagnostic use.

CAS	1204834-00-3
Purity	> 98%
Spacers	dPEG® Spacer is 28 atoms and 32.2 Å
Shipping	Ambient
Typical solubility properties (for additional information contact Customer Support)	Methylene chloride, Acetonitrile, DMAC, or DMSO.
Storage and handling	-20°C; Always let come to room temperature before opening; be careful to limit exposure to moisture and restore under an inert atmosphere; stock solutions can be prepared with dry solvent and kept for several days (freeze when not in use). dPEG® pegylation compounds are generally hygroscopic and should be treated as such. This will be less noticeable with liquids, but the solids will become tacky and difficult to manipulate, if care is not taken to minimize air exposure.

For research use only. Not intended for animal or human therapeutic or diagnostic use.