




## **IMPACT® AMEC RED SUBSTRATE, PEROXIDASE (HRP)**

**SKU:** SK-4285



---

### **DESCRIPTION**

ImmPACT AMEC Red HRP Substrate utilizes Vector Laboratories ImmPACT technology to generate a crisper, brighter red reaction product that is 5 -10 times more sensitive than conventional AEC substrates and can be used as a single label or as a second color for multiple antigen labeling.

#### **Features:**

- Optimized for IHC/ICC/ISH and blots
- Aqueous mounting
- Working solutions stable for up to 14 days at 4 °C
- Supplied as 120 ml of diluent and concentrated stock solutions in dropper bottles

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



## SPECIFICATIONS

<b>Chromogen Color</b>	Red
<b>Detection Enzyme</b>	Peroxidase
<b>Microscopy</b>	Brightfield, Spectral Imaging
<b>Mounting Medium Type</b>	Aqueous (Hardening), Aqueous (Non-Hardening)
<b>Unit Size</b>	120 ml
<b>Applications</b>	Immunohistochemistry / Immunocytochemistry, In situ hybridization, Blotting Applications, Elispot
<b>Additional Substrate Properties</b>	Contrast with Tissue Pigments, Multiple Labeling
<b>Safety Title</b>	WARNING: Cancer - <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a>
<b>Substrate Reaction Product</b>	Precipitating Reaction Product

## KIT COMPONENTS

- 120 ml ImmPACT AMEC Red Diluent
- 3 ml ImmPACT AMEC Red Reagent 1
- 3 ml ImmPACT AMEC Red Reagent 2

## TECHNICAL INFORMATION


ImmPACT AMEC Red chromogen can be used for both manual and automated staining methods. Stained slides must be aqueously mounted because the reaction product is soluble in organic solvents. Precipitated ImmPACT AMEC Red is stable for at least 2 years when mounted in VectaMount<sup>®</sup> AQ (H-5501).

ImmPACT AMEC Red consists of 120 ml of diluent and concentrated stock solutions of ImmPACT reagents in convenient dropper bottles. This substrate formulation can be used on nylon or PVDF membranes but is not recommended for use on nitrocellulose membranes.

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



### Relative Sensitivity of Substrates in Immunohistochemistry



\* Reaction product deposition not discreet or reliable for IHC. \*\* Longer incubation times increase sensitivity.

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



## Counterstain/Substrate Compatibility Table

This table is designed as a reference to determine the optimal counterstain/substrate combination for your application. Considerations should be given to tissue type, antigen unmasking protocol and other detection parameters to achieve the desired staining intensity.

Substrate	Catalog Number	Vector Hematoxylin and Hematoxylin QS H-3401 and H-3404	Vector Methyl Green H-3402	Vector Nuclear Fast Red H-3403
ImmPACT DAB (brown)	SK-4105	Excellent Contrast	Excellent Contrast	Fair Contrast
ImmPACT DAB EqV	SK-4103	Excellent Contrast	Excellent Contrast	Fair Contrast
DAB (brown)	SK-4100	Excellent Contrast	Excellent Contrast	Fair Contrast
DAB-Ni (gray-black)	SK-4100	Excellent Contrast	Fair Contrast *	Good Contrast
ImmPACT AEC (red)	SK-4205	Excellent Contrast	Counterstain Incompatibility **	Color Incompatibility
ImmPACT AMEC Red (red)	SK-4285	Excellent Contrast	Counterstain Incompatibility **	Color Incompatibility
AEC (red)	SK-4200	Excellent Contrast	Counterstain Incompatibility **	Color Incompatibility
TMB (blue)	SK-4400	Color Incompatibility	Counterstain Incompatibility	Excellent Contrast
ImmPACT VIP (purple)	SK-4605	Fair Contrast	Excellent Contrast	Poor Contrast
Vector VIP (purple)	SK-4600	Fair Contrast	Excellent Contrast	Poor Contrast
ImmPACT SG (blue-gray)	SK-4705	Poor Contrast	Good Contrast	Excellent Contrast
SG (blue-gray)	SK-4700	Poor Contrast	Good Contrast	Excellent Contrast
ImmPACT NovaRED (red)	SK-4805	Excellent Contrast	Excellent Contrast ***	Color Incompatibility
Vector NovaRED (red)	SK-4800	Excellent Contrast	Excellent Contrast ***	Color Incompatibility
ImmPACT Vector Red (magenta)	SK-5105	Excellent Contrast	Excellent Contrast	Color Incompatibility
Vector Red (magenta)	SK-5100	Excellent Contrast	Excellent Contrast	Color Incompatibility
Vector Black (black)	SK-5200	Excellent Contrast	Excellent Contrast *	Excellent Contrast
Vector Blue (blue)	SK-5300	Color Incompatibility	Good Contrast	Excellent Contrast
BCIP/NBT (indigo)	SK-5400	Color Incompatibility	Excellent Contrast *	Excellent Contrast

\* This substrate shows a slight decrease in sensitivity following the methyl green protocol. This decrease can be minimized by reducing the heat incubation and acetone rinse times in the methyl green protocol.

\*\* Substrate dissolves in acetone wash.

\*\*\* A slight color change in ImmPACT NovaRED™ and Vector NovaRED™ reaction product may be seen using methyl green.

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



### Enzyme Substrate Combinations

Recommended combinations of substrates and the recommended order in which they should be used.

Second Substrate		Alkaline Phosphatase			Peroxidase					
		ImmPACT Vector Red & Vector Red (magenta) SK-5105 & SK-5100	Vector Blue (blue) SK-5300	BCIP/NBT (indigo) SK-5400	ImmPACT VIP & Vector VIP (purple) SK-4605 & SK-4600	ImmPACT DAB, ImmPACT DAB EqV & DAB (brown) SK-4105, SK-4103, SK-4100	DAB-Ni (gray-black) SK-4100	ImmPACT NovaRED & Vector NovaRED (red) SK-4805 & SK-4800	ImmPACT SG & SG (blue-gray) SK-4705 & SK-4700	ImmPACT AEC, ImmPACT AMEC Red & AEC (red) SK-4205, SK-4285, SK-4200
First Substrate	Alkaline Phosphatase									
		ImmPACT Vector Red & Vector Red (magenta) SK-5105 & SK-5100		-	-	-	+	+	-	+
	Vector Blue (blue) SK-5300	+		-	+	+	+	+	+	+
	BCIP/NBT (indigo) SK-5400	+	-		+	+	+	+	+	+
Peroxidase	ImmPACT VIP & Vector VIP (purple) SK-4605, SK-4600	-	+	-		+	+	-	+	-
	ImmPACT DAB, ImmPACT DAB EqV & DAB (brown) SK-4105, SK-4103, SK-4100	+	+	+	+		-	-	+	+
	DAB-Ni (gray-black) SK-4100	+	-	-	+	+		+	-	-
	ImmPACT NovaRED & Vector NovaRED (red) SK-4805, SK-4800	-	+	+	-	+	+		+	-
	ImmPACT SG & SG (blue-gray) SK-4705, SK-4700	+	-	-	+	+	-	-		+
	ImmPACT AEC, ImmPACT AMEC Red & AEC (red) SK-4205, SK-4285, SK-4200	-	-	-	-	+	-	-	+	

+ Indicates good contrast      - Indicates incompatibility of substrates for various reasons

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz


## DOCUMENTS

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



- [Safety Data Sheet](#)
- [Enzyme Substrates Brochure](#)
- [Enzyme Substrate Combinations Chart](#)
- [Counterstain/Substrate Compatibility Chart](#)
- [User Guide](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**