




IMPRESS HRP HORSE ANTI-RABBIT IGG POLYMER DETECTION KIT, PEROXIDASE

SKU: MP-7401



DESCRIPTION

The ImmPRESS polymerized reporter enzyme staining system uses novel conjugation and micropolymer chemistries to create a highly sensitive, ready-to-use, one-step, non-biotin detection system for immunohistochemistry and immunocytochemistry staining. This unique micropolymer of highly active horseradish peroxidase (HRP) is attached to our affinity purified secondary antibodies, producing reagents with outstanding sensitivity and low background.

Features:

- Enhanced accessibility to nuclear and membrane antigens
- Ready-to-use in a convenient dropper bottle
- Shorter assay time
- Simplified multiple labeling
- 15 ml kit will stain approximately 75-150 sections and the 50 ml kit will stain approximately 250-500 sections

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Kit Contents:

- ImmPRESS (Peroxidase) Polymer Anti-Rabbit IgG Reagent (made in horse, ready-to-use)
- 2.5% Normal Horse Serum for blocking (ready-to-use)

SPECIFICATIONS

| | |
|-----------------------------|---------------------------------------------------------------------------|
| Detection Enzyme | Peroxidase |
| Format | Ready-to-Use |
| Isolation Technology | Micropolymer Reagents |
| Unit Size | 15 ml, 50 ml |
| Applications | Immunohistochemistry / Immunocytochemistry, ELISAs, In situ hybridization |
| Target Species | Rabbit |
| Conjugate | Micropolymer HRP |
| Source Species | Rabbit |


TECHNICAL INFORMATION

The peroxidase micropolymers of the ImmPRESS HRP polymer reagent limit steric interference and provide enhanced accessibility to the target, avoiding the disadvantages of other polymer systems that use large dextrans or other macromolecules as backbones. The result is crisp, strong staining of antibody targets, especially nuclear and membrane antigens (such as Ki67, estrogen receptor, bcl-2, CD3, CD8 and CD10) and greater sensitivity than other polymer systems.

Staining Procedure

The IHC immunostaining protocol is simple as shown in the diagram below. Following a blocking step with the diluted normal horse serum, sections are incubated with primary antibody. After a brief wash, the appropriate ImmPRESS Reagent is added to the sections and incubated for 30 minutes. Sections are again rinsed and the slides are developed with the peroxidase substrate of choice.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Consider Species Cross-Reactivity

When choosing the optimal detection system for your application, it is important to consider not only the species of the primary antibody but also the species of the tissue. If the species of the primary antibody and the species of the tissue are closely related (e.g. rat and mouse), the secondary antibody may bind to endogenous IgG in the tissue section leading to background. The following options minimize background staining in these instances:

- Use a secondary antibody specifically adsorbed to remove cross-reacting antibodies of closely-related species (e.g. ImmPRESS Anti-Mouse IgG, Rat Adsorbed).
- Use the M.O.M.[®] ImmPRESS Kit (Cat. No. MP-2400) for applications of mouse primary antibodies on mouse tissue.

Product review on Biocompare.com.

CITATIONS



Powered by Bioz © 2023 See more details on Bioz


DOCUMENTS

- [ImmPRESS Polymer & Automated Staining App. Note](#)
- [ImmPRESS Troubleshooting Guide](#)
- [User Guide](#)
- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.