


MB 660R NHS ESTER

SKU: FP-1661



Description

633/647



Laser
line

Cy5



Common
filter set

650



Excitation
max

665



Emission
max

MB™ 660R is a bright and photostable far-red dye that emits fluorescence at about 685 nm in the borderline spectral region between far-red and near-IR. Although the absorption maximum is at around 665 nm, this dye can be sufficiently excited by the 633 or 635 nm laser. MB™ 660R dye is water soluble and pH-insensitive from pH 4 to pH 10. MB 660R is a rhodamine-based dye, and like rhodamine dyes in general, it is exceptionally photostable (Figure 1, description tab). The superior photostability and excellent brightness of MB 660R make the dye an ideal choice for confocal microscopy and other demanding applications.


For research use only. Not intended for animal or human therapeutic or diagnostic use.

The NHS Ester reacts specifically and efficiently with a primary amine (e.g., side chain of lysine residues or aminosilane-coated surfaces) at pH 7-9 to form a stable, covalent amide bond. The NHS ester (or succinimidyl ester) is the most popular tool for conjugating dyes to the primary amines of protein or antibody (Lys), amine-modified oligonucleotides, and other amine-containing molecules.

MB™ 660R NHS Ester is an amine reactive, far-red emitting dye routinely used to label proteins or antibodies through the primary amines (Lys), amine-modified oligonucleotides, and other amine-containing biomolecules. The labeling occurs most efficiently at pH 7-9 and forms a stable, covalent amide bond.

MB™ 660R is a bright and photostable far-red dye that emits fluorescence at about 685 nm in the borderline spectral region between far-red and near-IR. Although the absorption maximum is at around 665 nm, this dye can be sufficiently excited by the 633 or 635 nm laser. MB™ 660R dye is water soluble and pH-insensitive from pH 4 to pH 10. MB 660R is a rhodamine-based dye, and like rhodamine dyes in general, it is exceptionally photostable (Figure 1). The superior photostability and excellent brightness of MB 660R make the dye an ideal choice for confocal microscopy and other demanding applications.


For research use only. Not intended for animal or human therapeutic or diagnostic use.



MB 660R dye spectrally is almost identical to [Alexa Fluor® 660](#) and [CF® 660R Dye](#) and can be used as a less expensive alternative to these dyes. Alexa Fluor® is a registered trademark of Thermo Fisher Scientific. CF® Dye is a registered trademark of Biotium.

Abs/Em Spectra

For research use only. Not intended for animal or human therapeutic or diagnostic use.



Specifications

Unit Size	1 mg, 5 mg, 25 mg, 100 mg
Reactivity	Primary amine
Abs/Em Maxima	665/690 nm
Extinction coefficient	92,000 cm ⁻¹ M ⁻¹
Solubility	Water, DMSO, DMF
Spectrally similar dyes	Alexa Fluor® 660, CF® 680R
Molecular weight	840.90 (protonated)
Storage Conditions	-20°C.
Shipping Conditions	Ambient temperature

For research use only. Not intended for animal or human therapeutic or diagnostic use.