



SUCCINYLATED WHEAT GERM AGGLUTININ (WGA), AGAROSE BOUND

SKU: AL-1023S



DESCRIPTION

Agarose bound, succinylated WGA is prepared using our affinity-purified lectins. This derivative has been reported to have properties distinct from the native lectin. Evidence suggests that Succinylated Wheat Germ agglutinin does not bind to sialic acid residues, unlike the native form, but retains its specificity toward *N*-acetylglucosamine. Using conjugates of the native lectin and the succinylated form can provide a system to distinguish between sialylated glycoconjugates and those containing only *N*-acetylglucosamine structures.

Features:

- Bead diameter ranges in size from 45-165 microns
- Matrix is stable in solutions at pH 3-11 as well as many organic solvents
- Immobilized lectins are prepared using affinity purified lectins
- Covalent attachment preserves lectin activity and minimizes conformational changes that might result in nonspecific or hydrophobic interactions

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



- Product supplied as a 1:1 suspension in buffer
- Inhibiting/Eluting Sugar: Chitin Hydrolysate; or 500 mM *N*-acetylglucosamine with salt and/or acid elution generally required; or Glycoprotein Eluting Solution (ES-5100)

SPECIFICATIONS

Molecular Weight	36
Extinction Coefficient	1.46
Formulation	10 mM HEPES, pH 7.5, 0.15 M NaCl, 20 mM GlcNAc, 0.08% sodium azide
Inhibiting or Eluting Sugar	Chitin Hydrolysate
Label Modifier Type	Lectins
Unit Size	2 ml
Storage Instructions	2-8 °C DO NOT FREEZE
Sugar Specificity	Terminal N-acetyl-containing glycans Wash gel thoroughly with buffer before use to remove sugar added to stabilize the lectin. Recommended product for eluting glycoconjugates bound to this agarose-lectin: Glycoprotein Eluting solution, Cat. No. ES-5100. Alternatively, 0.5 M N-Acetyl-D-Glucosamine (GlcNAc) can be used. For those glycoconjugates having very high affinity for WGA, it may be necessary to lower the pH of the eluting sugar solution to pH 3.0 with acetic acid and increase the concentration of GlcNAc. After use, wash the gel with several column volumes of buffered saline, then resuspend gel in buffered saline containing 0.08% sodium azide for storage.
Usage Summary	
Conjugate	Agarose

TECHNICAL INFORMATION

Agarose bound*, succinylated WGA is prepared using our affinity-purified lectins. Heat stable, cross-linked 4% agarose beads with a molecular weight exclusion limit of about 2×10^7 daltons are used as the solid-phase matrix to which the lectins are covalently coupled. The attachment

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



of the lectins to the beads is carefully controlled to preserve lectin activity and minimize conformational changes of the bound lectins that might result in nonspecific ionic or hydrophobic interactions. The technique we have developed to couple lectins to agarose beads inserts a hydrophilic spacer arm between the lectin and the matrix.

This coupling method provides several advantages over the traditional cyanogen bromide procedure:

- Maximum carbohydrate binding activity of the coupled lectins is retained
- Linkage is stable over a range of pH values
- Conjugated proteins are not leached off the beads by Tris or other routinely used buffers
- No residual charges are present after conjugation. This minimizes non-specific binding to the matrix

Our agarose bound lectins are supplied at a constant concentration of lectin per ml of settled beads. The concentration for each lectin is selected to achieve the highest glycoconjugate binding capacity per mg of lectin present in the beads. Each lot is tested for its binding capacity using glycoproteins known to bind the lectin. This provides a guideline for the user and assures the quality of our agarose bound lectins.

*1 mg lectin/ml gel

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Lectins in Histochemistry, ELISA, and Western Blot Applications](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES


For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Succinylated Wheat Germ Agglutinin \(WGA\), Agarose bound](https://vectorlabs.com/products/succinylated-wheat-germ-agglutinin-wga-agarose-bound/)

<https://vectorlabs.com/products/succinylated-wheat-germ-agglutinin-wga-agarose-bound/>