



# VECTASTAIN® ELITE ABC-HRP REAGENT, PEROXIDASE

**SKU:** PK-7100



---

## DESCRIPTION

The VECTASTAIN Elite ABC system is the most sensitive avidin/biotin-based peroxidase system and is approximately 5 times more sensitive than the original VECTASTAIN ABC Kit. It is supplied with the ready-to-use Elite ABC Kit reagent and does not include a biotinylated secondary antibody.

### Features:

- **Advanced avidin/biotin technology:** The Elite® ABC complex is smaller, very uniform, and highly active, allowing more accessibility for binding to a biotinylated target
- **Cost effective:** Higher sensitivity leads to lower cost per slide
- **Available in ready-to-use formats:** Prediluted, stabilized working solutions of Elite® ABC Kit reagents provide the same high sensitivity and low background as the traditional VECTASTAIN® Elite® ABC Kit reagents. Contains sufficient reagent to stain approximately 250-500 sections.

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



## SPECIFICATIONS

<b>Detection Enzyme</b>	Peroxidase
<b>Format</b>	Ready-to-Use
<b>Isolation Technology</b>	ABC Elite (Avidin/Biotin) Systems
<b>Unit Size</b>	50 ml
<b>Applications</b>	Immunohistochemistry / Immunocytochemistry, In situ hybridization, Blotting Applications, Elispot, ELISAs
<b>Target Species</b>	No antibody included
<b>Platform</b>	ABC Elite (Avidin/Biotin) Systems
<b>Source Species</b>	No antibody included

## KIT COMPONENTS

- 50 ml ready-to-use stabilized Elite® ABC reagent

## TECHNICAL INFORMATION

With more than 50,000 citations to its credit, the VECTASTAIN ABC kit remains widely popular. Based on the versatile avidin-biotin complex interaction, the system is modular, and along with our selection of secondary antibodies, can accommodate a wide array of primary antibody and tissue species. Our ABC kits are economical and continue to be a staple product in any immunohistochemistry (IHC) and immunocytochemistry (ICC) laboratory.

Peroxidase substrates produce sharp, dense precipitates with crisp localization. These characteristics, in conjunction with the high sensitivity and low background of the VECTASTAIN ABC systems, make the peroxidase enzyme a preferred choice in many applications. (eg. In neural tissue, the peroxidase system is often preferred because it gives more consistent labeling of both cell bodies and processes.)

### VECTASTAIN Elite ABC System

- **Advanced avidin-biotin technology:** The Elite ABC complex is smaller, very uniform, and highly active, allowing more accessibility for binding to a biotinylated target
- **Highest sensitivity, low background:** The VECTASTAIN Elite ABC system is the most

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



sensitive avidin-biotin complex-based peroxidase system. The Elite ABC series is approximately 5 times more sensitive than the original VECTASTAIN ABC Kit with the same low background.

- **Cost effective:** Higher sensitivity leads to lower cost per slide.
- Available without ( Standard kit) or with biotinylated species specific or universal secondary antibodies.
- **Available in Ready-To-Use (R.T.U.) formats:** Prediluted, stabilized working solutions of Elite ABC Kit reagents provide the same high sensitivity and low background as the traditional VECTASTAIN Elite ABC Staining Kit reagents.

## **VECTASTAIN Elite ABC Kit Component:**

The R.T.U. VECTASTAIN Elite ABC reagent consists of 50 ml stabilized, prediluted and preformed Elite ABC reagent. The R.T.U. kits are ideal for manual or automated tissue staining systems.

## **The Avidin Biotin Complex Method**

The VECTASTAIN ABC systems are extremely sensitive due to the form and number of active enzyme molecules associated with the preformed Avidin/Biotinylated enzyme Complex. This ABC complex takes advantage of two important properties of avidin: 1) an extraordinarily high affinity for biotin (over one million times higher than antibody for most antigens), and 2) four biotin-binding sites. These properties allow macromolecular complexes (ABCs) to be formed by mixing Avidin DH (Reagent A) with its paired biotinylated enzyme (Reagent B) prior to use. The ABC reagent, once formed, remains stable for many hours after formation and can be used for several days after preparation.

The VECTASTAIN ABC Reagent can be used to detect any molecule that is biotinylated. This property gives the avidin-biotin complex (ABC) method great versatility in the types of targets that can be detected as well as the types of applications in which it can be employed. Biotinylated primary antibodies, secondaries, lectins, neuronal tracers, nucleic acids, and ligands can be effectively visualized in applications such as:


- Tissue staining
- Multiple labeling (Multiplex IHC)
- Western blotting
- Southern and northern blotting
- In situ hybridization detection (ISH)
- Enzyme immunoassays (ELISA)
- Neuronal tracing

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



All applications benefit from the high sensitivity, low background, reproducibility, and economy of the VECTASTAIN ABC system.

### Using the ABC System:



## CHOOSING A VECTASTAIN ABC KIT

### 1) Choose an enzyme system:

The VECTASTAIN ABC Kits are available with a choice of two different enzyme systems

- Peroxidase (ABC, Elite® ABC, RTU or Kits)
- Alkaline Phosphatase (ABC-AP)

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



## 2) Choose the biotinylated secondary antibody:

To detect a biotinylated target, you will need only the ABC Reagent contained in the Standard Kit, followed by an appropriate substrate.

To detect an unlabeled primary antibody, you will need a biotinylated secondary antibody that binds to the primary antibody species chosen, the ABC Reagent, and an appropriate substrate. For example, to detect a primary antibody made in rabbit, the appropriate choice is a VECTASTAIN ABC Kit designated "Rabbit IgG".

For additional versatility and convenience, Vector Laboratories offers secondary antibody products that recognize more than one species of primary antibody. Our Universal Biotinylated Anti-Mouse/Rabbit IgG secondary antibody (Cat. No. BA-1400) is a cocktail of biotinylated anti-mouse IgG and anti-rabbit IgG, designed for use with both rabbit and mouse primary antibodies. This antibody should not be used to stain tissues from rodents or rabbits due to its reactivity with IgG present in the tissues of these species. The Universal Pan-Specific Anti-Mouse/Rabbit/Goat secondary antibody (Cat. No. BA-1300) (which is included in our VECTASTAIN Universal Kits) is a cocktail of biotinylated anti-mouse IgG, anti-rabbit IgG, and anti-goat IgG that recognizes mouse, rabbit, and goat primary antibodies. This secondary antibody should not be used to stain rat, mouse, or other rodent, rabbit, goat, bovine, or sheep tissue due to potential reactivity with endogenous IgG.

## 3) Choose a convenient format:

Finally, several kits are available in prediluted, stabilized formats for convenience and ease of use. See kits designated R.T.U. or Ready-To-Use.

## Consider Species Cross-Reactivity

When choosing the optimal detection system for your application, it is important to consider not only the species of the primary antibody but also the species of the tissue under examination. If the species of the primary antibody and the species of the tissue are closely related (e.g. rat and mouse), the biotinylated secondary antibody may cross react with endogenous IgG in the tissue section leading to background staining. The following three options minimize background staining in these instances:

- Biotinylate the primary antibody and detect with a Standard VECTASTAIN Elite ABC Kit, PK-6100 (only includes the avidin-biotin complex, or ABC reagent).
- Use a biotinylated secondary antibody specifically adsorbed to remove cross-reacting antibodies of closely-related species (e.g. biotinylated anti-mouse IgG, rat adsorbed, Cat. No. BA-2001).

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



- Use the M.O.M.<sup>®</sup> Immunodetection System for applications of mouse primary antibodies on mouse tissue.

## Substrates

After choosing the VECTASTAIN ABC Kit for your application, select a substrate that matches the enzyme system of the kit.

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [Avidin/Biotin Troubleshooting Guide](#)
- [User Guide](#)
- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**