




VECTASTAIN ELITE ABC UNIVERSAL PLUS KIT, PEROXIDASE (HORSE ANTI-MOUSE/RABBIT IGG)

SKU: PK-8200



DESCRIPTION

The VECTASTAIN Elite ABC system is the most sensitive avidin/biotin-based peroxidase system and is approximately 5 times more sensitive than the original VECTASTAIN ABC Kit with the same low background. Our Universal Biotinylated Anti-Mouse/Rabbit IgG secondary antibody included in this kit is a cocktail of biotinylated anti-mouse IgG and anti-rabbit IgG, designed for use with both rabbit and mouse primary antibodies. **This antibody should not be used to stain tissues from rodents or rabbits due to its reactivity with IgG present in the tissues of these species.**

Features:

- **Advanced avidin/biotin technology:** The Elite ABC complex is smaller, very uniform, and highly active, allowing more accessibility for binding to a biotinylated target
- **Cost effective:** Higher sensitivity leads to lower cost per slide

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



- **More complete detection systems:** Inclusion of BLOXALL[®] endogenous enzyme quenching solution and ImmPACT[®] DAB EqV substrate in the PLUS kit provides added convenience with optimized, matched reagents in the one product.

SPECIFICATIONS

Detection Enzyme	Peroxidase
Format	Ready-to-Use
Isolation Technology	ABC Elite (Avidin/Biotin) Systems
Unit Size	1 Kit
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization, Blotting Applications, Elispot, ELISAs
Target Species	Universal (Mouse/Rabbit)
Platform	ABC Elite (Avidin/Biotin) Systems

KIT COMPONENTS

- 50 ml ready-to-use stabilized Elite ABC reagent
- 50 ml ready-to-use universal, biotinylated anti-mouse/rabbit IgG
- 50 ml 2.5% normal horse serum
- 50 ml BLOXALL Endogenous Enzyme Blocking Solution
- ImmPACT DAB EqV Substrate (25 ml Reagent 1 + 25 ml Reagent 2)

TECHNICAL INFORMATION

With more than 50,000 citations to its credit, the VECTASTAIN ABC kit remains widely popular. Based on the versatile avidin-biotin complex interaction, the system is modular, and along with our selection of secondary antibodies, can accommodate a wide array of primary antibody and tissue species. Our ABC kits are economical and continue to be a staple product in any immunohistochemistry (IHC) and immunocytochemistry (ICC) laboratory.

Peroxidase substrates produce sharp, dense precipitates with crisp localization. These characteristics, in conjunction with the high sensitivity and low background of the VECTASTAIN ABC systems, make the peroxidase enzyme a preferred choice in many applications. (eg. In

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



neural tissue, the peroxidase system is often preferred because it gives more consistent labeling of both cell bodies and processes.)

VECTASTAIN Elite ABC Kit Components:

- 50 ml ready-to-use stabilized Elite® ABC reagent
- 50 ml ready-to-use universal, biotinylated anti-mouse/rabbit IgG
- 50 ml 2.5% normal horse serum
- 50 ml BLOXALL Endogenous Enzyme Blocking Solution
- ImmPACT DAB EqV Substrate (25 ml Reagent 1 + 25 ml Reagent 2)

Sufficient reagents to stain approximately 250-500 sections.

The Avidin Biotin Complex Method

The VECTASTAIN ABC systems are extremely sensitive due to the form and number of active enzyme molecules associated with the preformed Avidin/Biotinylated enzyme Complex. This ABC complex takes advantage of two important properties of avidin: 1) an extraordinarily high affinity for biotin (over one million times higher than antibody for most antigens), and 2) four biotin-binding sites. These properties allow macromolecular complexes (ABCs) to be formed by mixing Avidin DH (Reagent A) with its paired biotinylated enzyme (Reagent B) prior to use. The ABC reagent, once formed, remains stable for many hours after formation and can be used for several days after preparation.

Biotinylated primary antibodies, secondaries, lectins, neuronal tracers, nucleic acids, and ligands can be effectively visualized in applications such as:


- Tissue staining
- Multiple labeling (Multiplex IHC)
- Western blotting
- Southern and northern blotting
- In situ hybridization detection (ISH)
- Enzyme immunoassays (ELISA)
- Neuronal tracing

All applications benefit from the high sensitivity, low background, reproducibility, and economy of the VECTASTAIN ABC system.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Using the ABC System:



Choosing a VECTASTAIN ABC Kit

1) Choose an enzyme system:

The VECTASTAIN ABC Kits are available with a choice of two different enzyme systems

- Peroxidase (ABC, Elite[®] ABC, RTU or Kits)
- Alkaline Phosphatase (ABC-AP)

2) Choose the biotinylated secondary antibody:

To detect a biotinylated target, you will need only the ABC Reagent contained in the Standard Kit, followed by an appropriate substrate.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



To detect an unlabeled primary antibody, you will need a biotinylated secondary antibody that binds to the primary antibody species chosen, the ABC Reagent, and an appropriate substrate. For example, to detect a primary antibody made in rabbit, the appropriate choice is a VECTASTAIN ABC Kit designated “Rabbit IgG”.

For additional versatility and convenience, Vector Laboratories offers secondary antibody products that recognize more than one species of primary antibody. Our Universal Biotinylated Anti-Mouse/Rabbit IgG secondary antibody (Cat. No. BA-1400) is a cocktail of biotinylated anti-mouse IgG and anti-rabbit IgG, designed for use with both rabbit and mouse primary antibodies. This antibody should not be used to stain tissues from rodents or rabbits due to its reactivity with IgG present in the tissues of these species. The Universal Pan-Specific Anti-Mouse/Rabbit/Goat secondary antibody (Cat. No. BA-1300) (which is included in our VECTASTAIN Universal Kits) is a cocktail of biotinylated anti-mouse IgG, anti-rabbit IgG, and anti-goat IgG that recognizes mouse, rabbit, and goat primary antibodies. This secondary antibody should not be used to stain rat, mouse, or other rodent, rabbit, goat, bovine, or sheep tissue due to potential reactivity with endogenous IgG.

3) Choose a convenient format:

Finally, several kits are available in prediluted, stabilized formats for convenience and ease of use. See kits designated R.T.U. or Ready-To-Use.

Consider Species Cross-Reactivity

When choosing the optimal detection system for your application, it is important to consider not only the species of the primary antibody but also the species of the tissue under examination. If the species of the primary antibody and the species of the tissue are closely related (e.g. rat and mouse), the biotinylated secondary antibody may cross react with endogenous IgG in the tissue section leading to background staining. The following three options minimize background staining in these instances:

- Biotinylate the primary antibody and detect with a Standard VECTASTAIN Elite ABC Kit, PK-6100 (only includes the avidin-biotin complex, or ABC reagent).
- Use a biotinylated secondary antibody specifically adsorbed to remove cross-reacting antibodies of closely-related species (e.g. biotinylated anti-mouse IgG, rat adsorbed, Cat. No. BA-2001).
- Use the M.O.M.[®] Immunodetection System for applications of mouse primary antibodies on mouse tissue.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Substrates

After choosing the VECTASTAIN ABC Kit for your application, select a substrate that matches the enzyme system of the kit.

CITATIONS




Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Avidin/Biotin Troubleshooting Guide](#)
- [User Guide](#)
- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.