




## **Antigen Unmasking Solution, Citrate-Based**

H-3300-250

**Product Images** 



## **Short Description**

Vector Laboratories' Antigen Unmasking Solutions are highly effective at revealing antigens in formalin-fixed, paraffin embedded tissue sections when used in combination with a high temperature treatment procedure. This citrate-based solution (H-3300) is pH 6.0 and is supplied as convenient 100x concentrated stock, sufficient for preparation of 25 L of working solution.

## **Additional Information**

Unit Size	250 ml
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence

