




DAB Enhancing Solution

H-2200-30

[Product Images](#)



Short Description

DAB Enhancing Solution (H-2200) intensifies the DAB reaction product in stained sections. A ten second exposure to this solution darkens the reaction product.

Additional Information

Unit Size	30 ml
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization
Recommended Usage	*DAB Enhancing Solution is a ready-to-use reagent for toning DAB stained sections. Only visible DAB precipitate is toned, so little amplification is observed with the enhancing reagent.
Recommended Storage	2-8°C

