




Nuclear Fast Red Counterstain

H-3403-500

Product Images



Short Description

Nuclei stained with Vector® Nuclear Fast Red (H-3403) have a pink appearance. Vector Nuclear Fast Red is an excellent choice in multiple labeling applications when hematoxylin obscures the substrate colors.

Features:

- Good contrast with a variety of substrates
- Rapid staining procedure requiring only one step
- Suitable for use with non-aqueous mountants

Additional Information

Unit Size	500 ml
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization
Format	Ready-to-Use

