









VECTASHIELD® Antifade Mounting Medium

H-1000-10

Product Images





Short Description

VECTASHIELD Antifade Mounting Medium is a unique, stable formula for preserving fluorescence. This medium inhibits rapid photobleaching of fluorescent proteins and fluorescent dyes.

Features:

- Ideal refractive index
- Ready to use
- No warming necessary
- Can be stored without sealing for long term analysis
- Easy to use

The original VECTASHIELD Mounting Medium does not solidify, but remains a liquid on the slide and can be stored without sealing. If desired, coverslips can be sealed around the perimeter with nail polish or a plastic sealant. Mounted slides should be stored at 4 °C, protected from light.

Additional Information

Unit Size	10 ml
Applications	Immunofluorescence, In situ hybridization, Cellular Imaging, Super-Resolution Microscopy
Mounting	Aqueous (Non-Hardening)
Antifade	Yes
Counterstain	None

