






# VECTASHIELD® Antifade Mounting Medium with DAPI

## H-1200-10

[Product Images](#)



## Short Description

---

VECTASHIELD Antifade Mounting Medium with DAPI is a unique, stable formula for preserving fluorescence. This medium inhibits rapid photobleaching of fluorescent proteins and fluorescent dyes. VECTASHIELD Mounting Medium with DAPI does not solidify, but remains a liquid on the slide and can be stored without sealing.

### Features:

- Ideal refractive index
- Ready to use
- No warming necessary
- Can be stored without sealing for long term analysis
- Easy to use

## Additional Information

---

Unit Size	10 ml
Applications	Immunofluorescence, In situ hybridization, Cellular Imaging
Mounting	Aqueous (Non-Hardening)
Antifade	Yes
Counterstain	DAPI (blue fluorescence)
Cellular Stain	Nucleus

