





# VECTASHIELD® PLUS Antifade Mounting Medium

## H-1900

Product Images



## Short Description

---

VECTASHIELD PLUS protects against fading across the visible spectrum, even under far-red wavelengths. It is compatible with most commercially available fluorophores. This new formulation of non-setting media improves upon the original VECTASHIELD products. Specifically, VECTASHIELD PLUS exhibits no inherent background or toning and provides superior fluorophore signal retention across the spectrum, including far-red wavelengths.

### Features:

- No inherent toning or background
- Superior signal retention across the spectrum from blue to far red
- Mounted sections can be viewed in 1 hour
- Available in two convenient sizes
- Superior non-setting formulation


## Additional Information

---

Unit Size	2 ml, 10 ml
Applications	Immunofluorescence, In situ hybridization
Recommended Storage	2-8 °C
Mounting	Aqueous (Non-Hardening)
Antifade	Yes
Counterstain	None

### VECTASHIELD® PLUS Antifade Mounting Medium

VECTASHIELD PLUS Antifade Mounting Medium is a unique, stable formula for preserving fluorescence. VECTASHIELD PLUS Mounting Medium inhibits rapid photobleaching of fluorescent proteins and fluorescent dyes.




#### Features:

- Ideal refractive index
- Ready to use
- No warming necessary
- Can be stored without sealing for long term analysis
- Easy to use

The original VECTASHIELD PLUS Mounting Medium does not solidify, but remains a liquid on the slide and can be stored without sealing. If desired, coverslips can be sealed around the perimeter with nail polish or a plastic sealant. Mounted slides should be stored at 2-8 °C, protected from light.

### VECTASHIELD® PLUS Antifade Mounting Medium

VECTASHIELD PLUS Antifade Mounting Medium is a unique, stable formula for preserving fluorescence. VECTASHIELD PLUS Mounting Medium inhibits rapid photobleaching of fluorescent proteins and fluorescent dyes.



#### Features:

- Ideal refractive index
- Ready to use
- No warming necessary
- Can be stored without sealing for long term analysis
- Easy to use

The original VECTASHIELD PLUS Mounting Medium does not solidify, but remains a liquid on the slide and can be stored without sealing. If desired, coverslips can be sealed around the perimeter with nail polish or a plastic sealant. Mounted slides should be stored at 2-8 °C, protected from light.

