






VectaMount® AQ Aqueous Mounting Medium

H-5501-60

Product Images



Short Description

VectaMount AQ Mounting Medium (H-5501) is a hardsetting mounting medium developed for use with enzymatic substrates, such as AEC, ImmPACT® AEC, and ImmPACT AMEC Red, whose reaction products are soluble in alcohol or other organic solvents. In applications where aqueous mounting is preferred, VectaMount AQ is suitable for use with other substrates such as DAB, Vector® SG, BCIP/NBT, Vector Red, ImmPACT Vector Red, Vector Blue, ImmPACT DAB, ImmPACT DAB *EqV*, and ImmPACT SG.

Features:

- For aqueous mounting
- Compatible with many horseradish peroxidase and alkaline phosphatase substrates

Additional Information

Unit Size	60 ml
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization
Recommended Usage	VectaMount™ AQ Mounting Medium is a hard-setting aqueous-based mounting medium developed for use with enzymatic substrates such as AEC and ImmPACT™ AEC, whose reaction products are soluble in alcohol or other organic solvents. VectaMount AQ is also compatible with most peroxidase and alkaline phosphatase enzyme substrates, such as DAB, ImmPACT™ DAB, Vector® SG, ImmPACT™ SG, Vector® Red, Vector® Blue and BCIP/NBT.
Recommended Storage	Room Temperature*
Mounting	Aqueous (Hardening)
Antifade	No

