



TURBOLINK CATALYST BUFFER (10X)

SKU: S-2006-105



DESCRIPTION

TurboLINK Catalyst Buffer is used to catalyze the SoluLINK[®] bioconjugation reactions, which are based on aromatic hydrazone formation between HyNic-modified biomolecules and 4FB-modified biomolecules. TurboLINK greatly accelerates hydrazone formation — in the presence of TurboLINK Catalyst the conjugation reaction is typically complete in 1-2 hours. The active ingredient, aniline, is delivered at 100 mM concentration in Conjugation Buffer. A final concentration of 5-10 mM aniline in the conjugation reaction is sufficient for hydrazone catalysis. This reagent is particularly useful for conjugation of large biomolecules such as antibodies, oligonucleotides, proteins, enzymes, and surfaces.

SPECIFICATIONS

Unit Size	1.5 mL
Storage Instructions	2° - 8°C - Do Not Freeze

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Applications

In Situ Proximity Ligation

LEGAL NOTICE

SoluLINK® Bioconjugation For Research Use Only. Not for use in diagnostic procedures. For additional licensing restrictions, please see the license agreement at vectorlabs.com/solulink-research-license.

Products are for research use only, not for use in diagnostic or therapeutic procedures or for use in humans. Products are not for resale without express written permission of Seller. No license under any patent or other intellectual property right of Seller or its licensors is granted or implied by the purchase unless otherwise provided in writing.

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.