



AVIDIN DN

SKU: A-3100-1



DESCRIPTION

Avidin DN is isolated by an affinity chromatographic system that results in a form of avidin which, like Avidin DH, exhibits very low nonspecific binding to nucleic acids. Although the reasons for this property are not fully understood, this feature can be exploited when using biotin-avidin detection systems with biotin-labeled DNA or RNA probes for in situ or nitrocellulose hybridization techniques.

SPECIFICATIONS

Formulation	50 mM sodium bicarbonate, pH 8.2, preservative
Unit Size	1 mg
Storage Instructions	2-8 °C
Usage Summary	The recommended concentration range for use is 10-20 µg/ml. This product is designed for use with biotinylated nucleic acids.

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

Applications	In situ hybridization, Blotting Applications
Concentration	1.0 mg/ml
Conjugate	Unconjugated

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Avidin DN](#)

<https://vectorlabs.com/products/avidin-dn/>