




## **BLOXALL ENDOGENOUS BLOCKING SOLUTION, PEROXIDASE AND ALKALINE PHOSPHATASE**

**SKU:** SP-6000-100



---

### **DESCRIPTION**

BLOXALL Endogenous Peroxidase and Alkaline Phosphatase Blocking Solution is a dual enzyme blocking reagent that inactivates endogenous peroxidase, pseudoperoxidase, and alkaline phosphatase in formalin-fixed, paraffin-embedded tissue sections, frozen tissue sections, and cell preparations. This reagent is provided as a ready-to-use reagent in a convenient dropper bottle.

#### **Features: :**

- More effective than conventional blocking methods
- Single 10 minute incubation step
- 100 ml solution

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**



## SPECIFICATIONS

|                        |                                                                                                                                                                                                                                                                                                    |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Blocking Action</b> | Endogenous Peroxidase, Endogenous Alkaline Phosphatase                                                                                                                                                                                                                                             |
| <b>Unit Size</b>       | 100 ml                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>         | 2-8 °C                                                                                                                                                                                                                                                                                             |
| <b>Instructions</b>    |                                                                                                                                                                                                                                                                                                    |
| <b>Usage Summary</b>   | BLOXALL Endogenous Peroxidase and Alkaline Phosphatase Blocking Solution is a dual enzyme blocking reagent that inactivates endogenous peroxidase, pseudoperoxidase, and alkaline phosphatase in formalin-fixed, paraffin-embedded tissue sections, frozen tissue sections, and cell preparations. |
| <b>Applications</b>    | Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization                                                                                                                                                                                                              |

## TECHNICAL INFORMATION

Tissues can contain endogenous peroxidase, pseudoperoxidase, and/or alkaline phosphatase activity that will produce background staining if an alkaline phosphatase and/or peroxidase detection system and corresponding substrates are used.

## CITATIONS



Powered by Bioz © 2023 See more details on Bioz

## DOCUMENTS

- [User Guide](#)
- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**




[www.vectorlabs.com](http://www.vectorlabs.com)

Email: [customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

## GALLERY IMAGES



**For research use only. Not intended for therapeutic or diagnostic use in animals or humans.**

[BLOXALL Endogenous Blocking Solution, Peroxidase and Alkaline Phosphatase](https://vectorlabs.com/products/bloxall-endogenous-peroxidase-and-ap/)

<https://vectorlabs.com/products/bloxall-endogenous-peroxidase-and-ap/>