



CLEAVABLE BIOTIN ALKYNES

SKU: CCT-1454



DESCRIPTION

The sample pack contains 4 most often used alkyne-activated, cleavable biotin probes, Diazo, Dde, PC and DADPS Biotin Alkyne. These probes allow for efficient release of streptavidin-bound proteins under mild conditions. Even though all 4 cleavable moieties undergo efficient cleavage, the recovery of biopolymer or digests might vary from application to application requiring a preliminary testing of biotin probes to determine a best performing linker for a given application. Sample pack is a cost efficient option to profile all 4 cleavable linkers.

SPECIFICATIONS

Unit Size

Sample Pack

Storage Instructions

4C.

Shipping Instructions

Ambient temperature

DOCUMENTS

DATASHEETS:

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



- [Diazo Biotin Alkyne](#)
- [Dde Biotin Alkyne C](#)
- [PC Biotin Alkyne](#)
- [DADPS Biotin Alkyne](#)

SDS

- [Safety Data Sheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.