



CLEAVABLE BIOTIN AZIDES

SKU: CCT-1453



DESCRIPTION

This sample pack contains 4 of the most often used azide-activated, cleavable biotin probes: Diazo, Dde, PC, and DADPS Biotin Azide. These probes allow for the efficient release of streptavidin-bound proteins under mild conditions. Even though all 4 cleavable moieties undergo efficient cleavage, the recovery of biopolymer or digests might vary from application to application, requiring a preliminary testing of biotin probes to determine the best-performing linker for a given application. This sample pack is a cost-efficient option to profile all 4 cleavable linkers.

SPECIFICATIONS

Unit Size

Storage Instructions

Shipping Instructions

Sample Pack

4C.

Ambient temperature

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

DATASHEETS:

- [Diazo Biotin Azide](#)
- [Dde Biotin Azide](#)
- [PC Biotin Azide](#)
- [DADPS Biotin Azide](#)

SDS

- [Safety Data Sheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.