



CLICK-&-GO® MAGNETIC BEADS CONJUGATION KIT

SKU: CCT-1508



DESCRIPTION

Click-&-Go™ Magnetic Beads Conjugation Kit allows antibodies to be covalently attached to high quality 1 µm magnetic particles efficiently, quickly and easily using proprietary technology developed by scientists at Click Chemistry Tools. The conjugation reaction is initiated simply by the activation of the magnetic beads followed by briefly agitating the activated beads with an antibody which becomes attached (via lysine residues) to the magnetic particles. It takes only a few minutes to set up the conjugation, and the covalent conjugate is ready to use within 3 hours or less. Our kit provides enough material for 3 conjugations of 100 µg of antibody to 1.5 mg of the magnetic beads.

SPECIFICATIONS

Unit Size	1 Kit
Storage Instructions	4C
Number of Reactions	3
Shipping Instructions	Ambient temperature

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



DOCUMENTS

- [Safety Data Sheet](#)
- [User Guide](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.