

# CLICK-&-GO™ MAGNETIC BEADS CONJUGATION KITS

SKU: CCT-1508



## Description

Click-&-Go™ Magnetic Beads Conjugation Kit allows antibodies to be covalently attached to high quality 1  $\mu\text{m}$  magnetic particles efficiently, quickly and easily using proprietary technology developed by scientists at Click Chemistry Tools. The conjugation reaction is initiated simply by the activation of the magnetic beads followed by briefly agitating the activated beads with an antibody which becomes attached (via lysine residues) to the magnetic particles. It takes only a few minutes to set up the conjugation, and the covalent conjugate is ready to use within 3 hours or less. Our kit provides enough material for 3 conjugations of 100  $\mu\text{g}$  of antibody to 1.5 mg of the magnetic beads.

## Specifications

<b>Unit Size</b>	1 kit
<b>Number of Reactions</b>	3
<b>Storage Conditions</b>	4C
<b>Shipping Conditions</b>	Ambient temperature

---

**For research use only. Not intended for animal or human therapeutic or diagnostic use.**