



DAB ENHANCING SOLUTION

SKU: H-2200-30



DESCRIPTION

DAB Enhancing Solution (H-2200) intensifies the DAB reaction product in stained sections. A ten second exposure to this solution darkens the reaction product.

SPECIFICATIONS

Unit Size	30 ml
Storage Instructions	2-8°C
Usage Summary	*DAB Enhancing Solution is a ready-to-use reagent for toning DAB stained sections. Only visible DAB precipitate is toned, so little amplification is observed with the enhancing reagent.
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [User Guide](#)
- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.