



TEXAS RED™ AVIDIN D

SKU: A-2006-5



DESCRIPTION

Avidin D, Texas Red® labeled, is a highly fluorescent conjugate of Avidin D and sulforhodamine 101. Since the excitation and emission maxima are well separated from those of fluorescein, Texas Red® Avidin D can be employed with Fluorescein Avidin D or other fluorescein conjugates to simultaneously localize two antigens in the same tissue section.

SPECIFICATIONS

Color of Fluorescence Red

Format Concentrate

Formulation 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide

Maximum Emission 606-615 nm

Maximum Excitation 595-604 nm

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Unit Size	5 mg
Storage Instructions	2-8 °C
Usage Summary	For diluting this product, we recommend a HEPES- or bicarbonate-buffered saline solution, approximately pH 8.2. The recommended concentration range for use is 5-30 µg/ml.
Applications	Immunofluorescence, In situ hybridization, Flow Cytometry/Cell Separation
Concentration	2.5 mg active conjugate/mL
Conjugate	Texas Red

TECHNICAL INFORMATION

Texas Red[®] Avidin D is also ideal for flow cytometry using instruments equipped with dye lasers.

- Excitation maximum: 595 nm
- Emission maximum: 615 nm
- Color: Red

CITATIONS



Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



www.vectorlabs.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.

[Texas Red™ Avidin D](https://vectorlabs.com/products/texas-red-avidin-d/)

<https://vectorlabs.com/products/texas-red-avidin-d/>