




NUCLEAR FAST RED COUNTERSTAIN

SKU: H-3403-500



DESCRIPTION

Nuclei stained with Vector® Nuclear Fast Red (H-3403) have a pink appearance. Vector Nuclear Fast Red is an excellent choice in multiple labeling applications when hematoxylin obscures the substrate colors.

Features:

- Good contrast with a variety of substrates
- Rapid staining procedure requiring only one step
- Suitable for use with non-aqueous mountants

SPECIFICATIONS

Format Ready-to-Use

For research use only. Not intended for therapeutic or diagnostic use in animals or humans.



Unit Size 500 ml

Applications Immunohistochemistry / Immunocytochemistry, In situ hybridization

TECHNICAL INFORMATION

Optimize Counterstain Intensity

CITATIONS




Powered by Bioz © 2023 See more details on Bioz

DOCUMENTS

- [User Guide](#)
- [Counterstain/Substrate Compatibility Chart](#)
- [Safety Data Sheet](#)
- [Download CoA](#)
- [Datasheet](#)

GALLERY IMAGES



For research use only. Not intended for therapeutic or diagnostic use in animals or humans.