



---

## CERTIFICATE OF ANALYSIS

Product **BIOTINYLATED DEXTRAN AMINE (10Kd) TEXAS RED®**  
BDA-10KTR

Catalog No. SP-1140

Amount 10 mg

Lot No. ZB2111

Free Fluorochrome none

Approximate size 10,000 MW

Storage Conditions Refrigerate at 4 °C or, after reconstitution, store frozen.

Remarks:

Maximum Excitation @ 595-604 nm

Maximum Emission @ 606-615 nm

Reconstitute by addition of 0.4 ml of glass distilled water. The reconstituted product will have the following composition:

25 mM NaCl, pH 8.0. No preservative has been added.

<sup>®</sup>Texas Red is a Trademark of Molecular Probes