

VectaMount® Mounting Media & Substrate Compatibility Table

Substrate	Cat. No	VectaMount Express Mounting Medium (Non-Aqueous)	VectaMount Permanent Mounting Medium (Non-Aqueous)	VectaMount AQ Aqueous Mounting Medium
Horseradish Peroxidase (HRP)				
ImmPACT® DAB (brown)	SK-4105	✓	✓	✓
ImmPACT® DAB EqV (brown)	SK-4103	✓	✓	✓
DAB (brown)	SK-4100	✓	✓	✓
DAB-Ni (gray-black)	SK-4100	✓	✓	
ImmPACT® AEC (red)	SK-4205			✓
ImmPACT® AMEC Red (red)	SK-4285			✓
AEC (red)	SK-4200			✓
TMB (blue)	SK-4400		✓	
ImmPACT® VIP (purple)	SK-4605	✓	✓	
Vector® VIP (purple)	SK-4600	✓	✓	
ImmPACT® SG (blue-gray)	SK-4705	✓	✓	✓
Vector® SG (blue-gray)	SK-4700	✓	✓	✓
ImmPACT NovaRED® (red)	SK-4805	✓	✓	
Vector NovaRED® (red)	SK-4800	✓	✓	
Alkaline Phosphatase (AP)				
ImmPACT® Vector® Red (magenta)	SK-5105	✓	✓	✓
Vector® Red (magenta)	SK-5100	✓	✓	✓
Vector® Black (brown-black)	SK-5200		✓	
Vector® Blue (blue)	SK-5300	✓	✓	✓