




10x Casein Solution

SP-5020-250

Product Images



Short Description

10x Casein Solution is a general blocking agent for histochemical, nucleic acid and protein blotting, and other applications. It is supplied as 250 ml of a 10x concentrate.

Additional Information

Unit Size	250 ml
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs
Recommended Storage	2-8 °C
Recommended Usage	This product is designed to be used as a general blocking agent/diluent for nucleic acid and protein blotting applications. Prior to use, 10x Casein Solution should be diluted to 1x with distilled or deionized water. For prolonged storage of diluted Casein Solution, sodium azide should be added to give a final concentration of 0.08% (w/v) and the diluted solution refrigerated.
Blocking Action	Non-Specific Protein Blocking
Format	Concentrate

