



---

# Anti-DNP (Dinitrophenyl), Unconjugated

## SP-0603-1

[Product Images](#)

---



## Short Description

---

Dinitrophenyl (DNP) is not found endogenously in tissues and is an excellent alternative to biotin for use in applications requiring biotin-free systems. Probes or antibodies labeled with DNP can be detected using alkaline phosphatase labeled DNP or using the biotinylated antibody followed by a streptavidin or avidin system.

## Additional Information

---

Unit Size	1 mg
Applications	In situ hybridization, Blotting Applications
Concentration	1 mg active conjugate/ml
Recommended Storage	2-8 °C; Store frozen for long term storage
Solution	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide.
Detection Target	Dinitrophenyl (DNP)
Conjugate	Unconjugated
Host Species	Rabbit
Format	Concentrate

