





# Antigen Unmasking Solution, Tris-Based

## H-3301-250

Product Images



## Short Description

---

Vector Laboratories' Antigen Unmasking Solutions are highly effective at revealing antigens in formalin-fixed, paraffin embedded tissue sections when used in combination with a high temperature treatment procedure. This Tris-based solution (H-3301) is pH 9.0 and is supplied as a convenient 100x concentrated stocks sufficient for preparation of 25 L of working solution.

## Additional Information

---

Unit Size	250 ml
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence

