




# BLOXALL® Endogenous Blocking Solution, Peroxidase and Alkaline Phosphatase

## SP-6000-100

### Product Images



## Short Description

---

BLOXALL Endogenous Peroxidase and Alkaline Phosphatase Blocking Solution is a dual enzyme blocking reagent that inactivates endogenous peroxidase, pseudoperoxidase, and alkaline phosphatase in formalin-fixed, paraffin-embedded tissue sections, frozen tissue sections, and cell preparations. This reagent is provided as a ready-to-use reagent in a convenient dropper bottle.

### Features:

- More effective than conventional blocking methods
- Single 10 minute incubation step
- 100 ml solution

## Additional Information

---

Unit Size	100 ml
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Neurobiology
Recommended Storage	2-8 °C
Recommended Usage	BLOXALL Endogenous Peroxidase and Alkaline Phosphatase Blocking Solution is a dual enzyme blocking reagent that inactivates endogenous peroxidase, pseudoperoxidase, and alkaline phosphatase in formalin-fixed, paraffin-embedded tissue sections, frozen tissue sections, and cell preparations.
Blocking Action	Endogenous Peroxidase, Endogenous Alkaline Phosphatase

