



---

# Bovine Serum Albumin (BSA), Biotinylated

## B-2007-10

[Product Images](#)

---



## Short Description

---

This product is produced from our ultrapure immunohistochemical grade bovine serum albumin. This biotinylated protein can be used for measuring binding capacity of avidin or streptavidin matrices, as a reagent for amplifying biotin-avidin detection systems, or as a positive control in dot blot assays.

## Additional Information

---

Unit Size	10 mg
Applications	Blotting Applications, Affinity Chromatography
Solution	10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% azide.
Recommended Storage	2-8 °C
Concentration	5 mg/ml
Conjugate	Biotinylated

