



# Bovine Serum Albumin (BSA), Immunohistochemical Grade

## SP-5050-500

[Product Images](#)



## Short Description

---

This ultrapure grade of bovine serum albumin (BSA) can be used as a diluent or a blocking agent. It is free of impurities present in other grades of BSA which can introduce artifacts or increase background in immunohistochemical staining, ELISAs, or blot development.

## Additional Information

---

Unit Size	500 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs, Flow Cytometry/Cell Separation
Recommended Storage	2-8 °C
Blocking Action	Non-Specific Protein Blocking
Format	Lyophilized

