



Goat Anti-Human Antibody, Kappa Chain, Unconjugated

AI-3060-.5

[Product Images](#)



Short Description

Unconjugated Goat Anti-Human Kappa Chain Antibody binds to human primary antibodies for immunohistochemistry and many other applications. The unlabeled form can be used as a capture antibody in enzyme immunoassays or radioactively labeled for radioimmunoassays. This kappa chain specific antibody has virtually no cross-reactivity with other immunoglobulin classes or other heavy or light chains.

Additional Information

Unit Size	0.5 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs
Concentration	1 mg/ml
Recommended Storage	2-8 °C; Store frozen for long term storage
Solution	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose
Target Species	Human
Conjugate	Unconjugated
Host Species	Goat
Format	Concentrate

