



Goat Anti-Human IgG Antibody, Gamma Chain Specific, Biotinylated

BA-3080-.5

[Product Images](#)



Short Description

Biotinylated Goat Anti-Human IgG Antibody, gamma chain specific, is used in an avidin-biotin or streptavidin-biotin detection system. These antibodies can be used for tissue and cell staining, ELISAs, and blots. This gamma chain specific antibody has virtually no cross-reactivity with other immunoglobulin classes or other heavy or light chains.

Additional Information

Unit Size	0.5 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs
Concentration	0.5 mg active conjugate/ml
Recommended Storage	2-8 °C; Store frozen for long term storage
Solution	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 1 mg/ml bovine serum albumin.
Recommended Usage	The recommended concentration range for use is 2-10 µg/ml. If this biotinylated antibody is to be used in tissues, which may contain cross-reacting endogenous immunoglobulins, dilution of this biotinylated antibody may be made in buffers containing 2% normal serum from the same species as the tissue.
Target Species	Human
Conjugate	Biotinylated
Host Species	Goat
Format	Concentrate

