



---

# Goat Anti-Human IgG Antibody, Gamma Chain, Unconjugated

## AI-3080-.5

[Product Images](#)

---



## Short Description

---

Unconjugated Goat Anti-Human IgG, gamma chain specific, binds to human primary antibodies for immunohistochemistry and many other applications. This gamma chain specific antibody has virtually no cross-reactivity with other immunoglobulin classes or other heavy or light chains.

### Features:

- Affinity purified, ultrapure, high affinity antibody
- Thoroughly adsorbed against serum and immunoglobulins from potentially interfering species
- Ready for iodination, fluorochrome labeling, or enzyme conjugations
- Can be employed as a capture antibody in enzyme immunoassays or in other assays requiring carrier-free immunoglobulins

## Additional Information

---

Unit Size	0.5 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs
Concentration	1 mg/ml
Recommended Storage	2-8 °C; Store frozen for long term storage
Solution	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose
Target Species	Human
Conjugate	Unconjugated
Host Species	Goat
Format	Concentrate

