



---

# Goat Anti-Human Kappa Chain Antibody, Biotinylated

## BA-3060-.5

[Product Images](#)

---



## Short Description

---

Biotinylated Goat Anti-Human Kappa Chain Antibody is a chain-specific antibody that can be used for tissue and cell staining, ELISAs, and blots. This kappa chain specific antibody has virtually no cross-reactivity with other immunoglobulin classes or other heavy or light chains.

### Features:

- Biotinylated to ensure the maximum degree of labeling without compromising the specificity or affinity of the antibody

## Additional Information

---

Unit Size	0.5 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs
Recommended Storage	2–8 °C
Recommended Usage	Reconstitute by adding 1.0 ml water. The resulting solution will have the following composition: 10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 1 mg/ml bovine serum albumin. The recommended concentration range for use is 2-10 µg/ml. If this biotinylated antibody is to be used in tissues, which may contain cross-reacting endogenous immunoglobulins, dilution of this biotinylated antibody may be made in buffers containing 2% normal serum from the same species as the tissue.
Target Species	Human
Conjugate	Biotinylated
Host Species	Goat
Format	Concentrate

