



Goat Anti-Rabbit IgG Antibody (H+L), Texas Red™

TI-1000-1.5

[Product Images](#)



Short Description

The Texas Red Goat Anti-Rabbit IgG Antibody can be used for immunofluorescence and other applications. Optimal F/P ratios have been established for each conjugate to ensure maximum fluorescence with minimal background staining.

Features:

- Recognizes both heavy and light chains (H+L)
- Optimally labeled with Texas Red to provide the brightest label for fluorescence microscopy
- Supplied in solution
- Excitation: 595 nm
- Emission: 615 nm
- Color: Red

Additional Information

Unit Size	1.5 mg
Applications	Immunofluorescence, In situ hybridization, Blotting Applications, Flow Cytometry/Cell Separation
Concentration	1.5 mg active conjugate/ml
Recommended Storage	2-8 °C
Solution	10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.
Maximum Emission	606-615 nm
Maximum Excitation	595-604 nm
Recommended Usage	The recommended concentration range for use is 5-20 µg/ml. If this antibody is to be used in tissues which may contain cross-reacting endogenous immunoglobulins, dilution of this antibody may be made in buffers containing 2 % normal serum from the same species as the tissue
Target Species	Rabbit
Conjugate	Texas Red
Color of Fluorescence	Red
Host Species	Goat
Format	Concentrate

