





# Hematoxylin Counterstain

## H-3401-500

### Product Images



## Short Description

---

Vector Hematoxylin is based on Gill's formulation, it produces blue stained nuclei with crisp detail.

### Features:

- Alcohol and mercury free
- Progressive stain formula
- Provides excellent color contrast with most commonly used peroxidase and alkaline phosphatase substrates
- Suitable for use with non-aqueous and aqueous mountants

## Additional Information

---

Unit Size	500 ml
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization
Format	Ready-to-Use

