









Hematoxylin QS Counterstain

H-3404-100

Product Images



Short Description

Vector Hematoxylin QS (H-3404), a modification of Mayer's hematoxylin developed especially for immunocytochemistry, provides crisp blue nuclear staining. This formula is ready-to-use without filtration, requires no separate bluing step, and stains in less than 45 seconds.

Features:

- Provides excellent color contrast with most commonly used peroxidase and alkaline phosphatase substrates
- Suitable for use with nonaqueous and aqueous mountants
- Mercury-free

Additional Information

Unit Size	100 ml
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization
Format	Ready-to-Use

