



Horse Anti-Mouse IgG Antibody (H+L), Peroxidase

PI-2000-1

[Product Images](#)



Short Description

The HRP Horse Anti-Mouse IgG Antibody binds to mouse primary antibodies for tissue staining, ELISAs, and blots. It can be detected using a horseradish peroxidase substrate.

Features:

- Recognizes both heavy and light chains (H+L)
- Cross-linked with the highest specific activity horseradish peroxidase
- Supplied in solution

Additional Information

Unit Size	1 mg
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization, Blotting Applications, Elispot, ELISAs
Concentration	1 mg active conjugate/ml
Recommended Storage	2-8 °C, DO NOT FREEZE
Solution	10 mM HEPES, 0.15 M NaCl, pH 7.5, preservative
Recommended Usage	Recommended diluent for microtitration plate assays: 10 mM phosphate, 0.15 M NaCl, pH 7.5, 0.1% Tween 20. Working concentration of HRP-anti mouse IgG is 0.5-5.0 µg/ml.
Target Species	Mouse
Conjugate	Peroxidase
Host Species	Horse
Format	Concentrate
Safety Title	Cancer

