



# Horse Anti-Mouse IgG Antibody (H+L), Unconjugated

## AI-2000

[Product Images](#)



## Short Description

---

Unconjugated Horse Anti-Mouse IgG Antibody binds to mouse primary antibodies for immunohistochemistry and many other applications.

### Features:

- Recognizes both heavy and light chains (H+L)
- Ready for iodination, fluorochrome labeling, or enzyme conjugations
- Can be employed as a capture antibody in enzyme immunoassays or in other assays requiring carrier-free immunoglobulins
- Supplied in liquid format

## Additional Information

---

Unit Size	1.5 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs
Concentration	3 mg/ml
Recommended Storage	2-8 °C; Store frozen for long term storage
Solution	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose
Target Species	Mouse
Conjugate	Unconjugated
Host Species	Horse
Format	Concentrate

