



ImmPrint™ Permanent Marking Pen

H-6100

Product Images



Short Description

The ImmPrint Permanent Marking Pen is a solvent-resistant permanent marking pen that is designed for writing on frosted glass microscope slides, tissue cassettes, and most hard surfaces.

Features:

- High density, fast-drying, black ink is resistant to most organic solvents encountered in histological applications such as alcohols and clearing agents
- Has a smooth writing tip that resists drying out

Additional Information

Unit Size	5 pen set
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization

