





# Methyl Green Counterstain

## H-3402-500

### Product Images



# Short Description

Vector Methyl Green stains cell nuclei light green. It is an excellent alternative in multiple labeling when hematoxylin obscures the substrate colors.

**Features:**

- Simple, two-step procedure
- Suitable for use with a wide range of enzyme substrates
- Suitable for use with non-aqueous mountants

# Additional Information

Unit Size	500 ml
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization
Format	Ready-to-Use

