



# Normal Goat Serum Blocking Solution, 2.5%

## S-1012-50

[Product Images](#)



## Short Description

---

2.5% Normal Goat Serum can be used for blocking non-specific binding or as an antibody diluent.

### Features:

- Pooled samples collected from healthy adult animals
- Heat-treated and centrifuged to remove precipitates and then filtered
- Tested with appropriate antibody to check for possible cross-reactivities
- Ready-to use 2.5% solution (v/v) in buffer containing 0.08% sodium azide
- Same serum preparation that is included in the ImmPRESS® Detection Systems

## Additional Information

---

Unit Size	50 ml
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence
Blocking Action	Non-Specific Protein Blocking
Format	Ready-to-Use

