



Rabbit Anti-Goat IgG Antibody (H+L), Unconjugated

AI-5000-1.5

[Product Images](#)



Short Description

Unconjugated Rabbit Anti-Goat IgG Antibody binds to goat primary antibodies for immunohistochemistry and many other applications. It is also suitable for use with sheep and bovine IgG primary antibodies

Features:

- Ultrapure, high affinity antibody
- Thoroughly adsorbed against serum and immunoglobulins from potentially interfering species
- Recognizes both heavy and light chains (H+L)
- Ready for iodination, fluorochrome labeling, or enzyme conjugations
- Can be employed as a capture antibody in enzyme immunoassays or in other assays requiring carrier-free immunoglobulins

Additional Information

Unit Size	1.5 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization, Blotting Applications, Elispot, ELISAs
Concentration	3 mg/ml
Recommended Storage	2-8 °C; Store frozen for long term storage
Solution	10 mM sodium phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide, 20 mg/ml sucrose
Target Species	Goat
Conjugate	Unconjugated
Host Species	Rabbit
Format	Concentrate

