



Rabbit IgG, Control Antibody

I-1000-5

[Product Images](#)



Short Description

This IgG preparation is intended for use as a control for primary antibodies made in rabbit. This antibody has been purified from pooled serum of healthy adult animals and contains a spectrum of the IgG subclasses present in serum.

Additional Information

Unit Size	5 mg
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, Blotting Applications
Recommended Storage	2-8 °C; Store frozen for long term storage
Recommended Usage	Reconstitute by adding 1 ml water. The resulting solution will have the following composition: 10 mM phosphate, 0.15 M NaCl, pH 7.8, 0.08% sodium azide, 20 mg/ml sucrose. ALLOW THE PROTEIN TO DISSOLVE WITHOUT VIGOROUS AGITATION.
Conjugate	Unconjugated
Host Species	Rabbit

