




# S-4FB Linker (DMF Soluble)

## S-1004

### Product Images



## Short Description

---

The S-4FB heterobifunctional crosslinker is fundamental to the SoluLINK bioconjugation technology. S-4FB reacts with primary amines on proteins (through lysine amino acids) or amino-modified oligonucleotides or surfaces, introducing a 4FB linker that forms stable covalent conjugates with biomolecules possessing HyNic linkers.

## Additional Information

---

Label/Modifier Type	4FB
Reactivity	HyNic
Recommended Storage	Desiccated: -15° to -25°C
Applications	Antibody Labeling, Antisense/RNAi, Aptamers
Other Name(s)	SoluLINK Bioconjugation

## Products in this set

---



**S-4FB Linker (DMF Soluble)**



**S-4FB Linker (DMF Soluble)**

