




Sulfo ChromaLINK® Biotin (Water Soluble)

B-1007

[Product Images](#)



Short Description

ChromaLINK Biotin contains a UV-traceable chromophore based on ChromaLINK technology to enable reproducibility in your biotinylation process. Now you can measure the degree of biotinylation in minutes, not hours, without the standard curves required for HABA/avidin and fluoro-reporter assays. With a simple, non-destructive UV scan, you can quantify biotin incorporation and ensure reproducible production of consistent batches. This biotinylation reagent is water soluble, which means no DMF or other organic solvent is required to dissolve the compound prior to labeling.

Additional Information

Label/Modifier Type	biotin
Reactivity	Streptavidin
Recommended Storage	Desiccated: -15° to -25°C
Applications	Antibody Labeling, Aptamers
Other Name(s)	SoluLINK Bioconjugation

Products in this set



Sulfo ChromaLINK® Biotin (Water Soluble)



Sulfo ChromaLINK® Biotin (Water Soluble)

