



VECTASTAIN® Elite® ABC-HRP Reagent, Peroxidase, R.T.U.

PK-7100

[Product Images](#)



Short Description

The VECTASTAIN Elite ABC system is the most sensitive avidin/biotin-based peroxidase system and is approximately 5 times more sensitive than the original VECTASTAIN ABC Kit. It is supplied with the ready-to-use Elite ABC Kit reagent and does not include a biotinylated secondary antibody.

Features:

- **Advanced avidin/biotin technology:** The Elite® ABC complex is smaller, very uniform, and highly active, allowing more accessibility for binding to a biotinylated target
- **Cost effective:** Higher sensitivity leads to lower cost per slide
- **Available in ready-to-use formats:** Prediluted, stabilized working solutions of Elite® ABC Kit reagents provide the same high sensitivity and low background as the traditional VECTASTAIN® Elite® ABC Kit reagents.

Contains sufficient reagent to stain approximately 250-500 sections.

Additional Information

Unit Size	50 ml
Applications	Immunohistochemistry / Immunocytochemistry, In situ hybridization, Blotting Applications, Elispot, ELISAs
Technology	ABC Elite (Avidin/Biotin) Systems
Detection Enzyme(s)	Peroxidase
Target Species	No antibody included
Format	Ready-to-Use

